

GLODERM

SKIN BOOSTERS | SKIN FILLERS

NATURAL | SMOOTH | SAFE



GLODERM I for natural beauty

GLODERM is a new and innovative range of smooth hyaluronic acid gels for enhancing natural beauty

GLODERM RANGE | FOR FACIAL REJUVENATION AND LIFTING



Skin Boosters
best suited for skin rejuvenation

Skin Fillers
tailored to smooth wrinkles
restore facial contours and volume

GLODERM 20 | GLODERM 30

- formulated with natural, non-crosslinked hyaluronic acid
- tailored for rejuvenation of mid-age, mature and photo-damaged skin
- versatile tools to improve skin hydration, structure and condition



20

HA: 20 mg/ml
natural, non crosslinked

- rejuvenation of global face, neck, decollete and hands
- mid-age, dry and thin skin
- lips hydration, smoothing of fine lines



30

HA: 30 mg/ml
natural, non crosslinked

- rejuvenation of global face
- mature and photo-damaged skin
- improvement of acne scars

GLODERM 20L | GLODERM 30L

- monophasic, smooth hyaluronic acid gels
- low crosslinked
- fully hydrated
- lifting of skin depressions
- restoring of facial contours and fullness
- immediate aesthetic results that last



20L

HA: 20 mg/ml
natural, crosslinked

- filling of moderate wrinkles
- lips reshaping



30L

HA: 30 mg/ml
natural, crosslinked

- filling of deep wrinkles and folds
- restoring of facial contours and volume

GLODERM FILLERS I

The science behind

GLODERM Fillers:


- formulated with Ion-Bonding Self Crosslinking Technology (IBS-pat.pend.)
- safe and reproducible
- smooth and easy to inject



GLODERM FILLERS | NATURAL AND SAFE



non cross-linked HA



cross-linked HA

- Network of HA chains with high density of natural cross-linking (flexible entanglement) and low density of synthetic cross-links

GLODERM FILLERS | RANGE OF INDICATIONS

Glabella
Gloderm 20L

Nasolabial folds
Gloderm 20L/30L

Oral commissures
Gloderm 30L

Jawline augmentation
Gloderm 30L

Tear trough
Gloderm 20L

Cheek augmentation
Gloderm 30L

Marionette lines
Gloderm 20L

Lip enhancement
Gloderm 20L

Chin augmentation
Gloderm 30L

GLODERM RANGE | NATURAL AND SAFE BEAUTY ENHANCERS

	GLODERM 20	GLODERM 30	GLODERM 20L	GLODERM 30L
HA Concentration (mg/ml)	20	30	20	30
Level of cross-linking	○○○	○○○	●●○	●●●
Viscosity	Intermediate	High	Intermediate	High
Suggested injection depth	Superficial dermis	Superficial to mid dermis	Mid dermis	Deep dermis
Lifting capacity	---	---	Moderate	High
Pack	1 x 1 ml	1 x 1 ml	1 x 1 ml	1 x 1 ml
Suggested Indications	<ul style="list-style-type: none"> · Rejuvenation of global face, neck, décolletage, hands · Lips hydration · Smoothing of fine 	<ul style="list-style-type: none"> · Rejuvenation of global face · Improvement of acne scars 	<ul style="list-style-type: none"> · Nasolabial folds · Glabella · Marionette lines · Lip contour and volume 	<ul style="list-style-type: none"> · Deep nasolabial folds · Oral commissures · Cheek, jawline and chin augmentation
Needle size	30 G	30 G	27 G	27G



Marketed By:

Biotech Healthcare Holding GmbH
Obergrundstrasse 17, 6002 Luern
Switzerland